
Patient interface selection



Ventilator compatibility											
Single limb circuit	EE elbows (leak 1, 2)	•	EE elbow	EE elbow	EE elbow	EE elbow	•	•	•	•	
Dual limb circuit	SE elbow		SE elbow	SE elbow	SE elbow	SE elbow	•	•			
Patient treatments											
Bronchoscopy procedure	Bronchoscopy elbow		Bronchoscopy elbow								
Aerosol medication delivery			NIVO elbows								
Patient characteristics											
Claustrophobia	•	•								•	
Facial irregularities	•	•								•	
With no teeth	•	•	•	•	•	•	•	•	•	•	
Easier mouth access			with CapStrap	•	with CapStrap		•	•	•	•	
Mouth breathing	•	•	•	•	•	•					
Application indications											I
Rapid application	•	•	with CapStrap	•	with CapStrap			•			
Frequent NIV patient	•		•	•	•	•	•	•		•	
Post-surgical/recovery			•	•	•	•					
Headgear type											I
Four-point	•	•	•		•		•				
Four-point with chin strap										•	
CapStrap			•	•	•			•			
Three-point						•					

Elbows



EE, Leak 1
(Respironics BIPAP devices)



EE, Leak 1 shown
with NIVO nebulizer
(Respironics BIPAP devices)



EE, Leak 2
(Respironics BIPAP devices)



(ICU ventilators)



SE shown
with NIVO nebulizer
(ICU ventilators)



Bronchoscopy elbow
(Respironics BIPAP &
ICU ventilators)

© 2011 Koninklijke Philips Electronics N.V. All rights are reserved. Respironics Total, PerforMax, PerformaTrak, Image3, Contour Deluxe, and CapStrap are trademarks of Respironics, Inc. and its affiliates.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

For more information contact your Philips representative.



Bonnet